

Orthopaedic

Product Range Guide



Edition 3

The conservative management of soft tissue injuries, fractures and orthopaedic conditions requires a broad range of high-quality, clinically effective products that are easy to use, cost-effective and maintain patient quality of life.

This range provides an all-encompassing portfolio to meet the treatment needs of patients requiring primary and/or secondary immobilisation through to post-healing rehabilitation.



Casting and Splinting Workshops

Essity has an ongoing commitment to educational support for clinicians involved in the provision of casting and splinting services in the healthcare environment. Training is offered in the form of in-services or educational workshops (with hands-on activity) to ensure all clinicians understand the appropriate and cost effective use of our products.

There is also a range of educational resources including application posters, instructional videos, application manuals and usage guides to assist in treatment management.


If you would like more detailed information about the education services on offer, please contact Essity Customer Service on 1300 276 633.

Contents

	Non-invasive fracture treatment pathway	4		Cast Protection	
				Cast Shoes	27
	Primary Immobilisation			Orthopaedic Soft Goods	
	Dynacast® Prelude	5		Actimove® Post-Op ROM Knee Brace	28
	Gypsona®	6		Actimove® Genu Knee Immobiliser	29
	Gypsona® Extra 3	7		Actimove® Genu Tri-Panel Knee Immobiliser	30
	Specialist®	8		Actimove® Lower Leg Walker	31
	Tensoplast® Skin Traction Kit	9		Actimove® Solea	32
	Secondary Immobilisation			Actimove® Shoulder Abduction Sling	33
	Delta-Lite® Plus	10		Actimove® Umerus Comfort	34
	Delta-Cast® Elite	11		Actimove® Umerus	35
	Delta-Cast® Prints	12		Actimove® Mitella Comfort	36
	Delta-Cast® Soft	13		Actimove® Mitella	37
	Undercast Materials			Actimove® Mitella for Kids	38
	Soffban® Natural	14		Actimove® Sling	39
	Soffban® Synthetic	15		Actimove® Sling for Kids	40
	Artiflex® Soft	16		Actimove® TaloMotion	41
	Velband™	17		Actimove® GenuMotion	42
	AquaCast® Liner	18		Actimove® ManuMotion	43
	Delta-Dry®	19		Actimove® Rhizo Forte	44
	Protouch® Stockinette	20		Actimove® Manus Forte	45
	Delta Terry-Net™ S	21		Actimove® Manus Forte Plus	46
	Cast Removal Equipment			Actimove® EpiMotion	47
	deSoutter Cast Saw CC4	22		Actimove® Cervical	48
	deSoutter Cast Saw CC5	23		Actimove® Vertebrace	49
	deSoutter Battery Powered Saw – CSB-100	24			
	deSoutter Battery/Mains Powered Saw CSP-201	25			
	Cast Room Accessories	26			

Please Note: Not all products available in New Zealand

Non-invasive fracture treatment pathway

Treatment Pathway	Objective	Considerations	Product Options
 <p>Primary Immobilisation</p>	<ul style="list-style-type: none"> Reduce and stabilise fracture Manage soft tissue injury Reduce oedema 	<ul style="list-style-type: none"> Pain Oedema Fracture stability Wound management 	<ul style="list-style-type: none"> Gypsona® Gypsona® Extra 3 Specialist® Dynacast® Prelude Tensoplast® Traction Kits Actimove® Supports Elastolite®
 <p>Secondary Immobilisation</p>	<ul style="list-style-type: none"> Maintain alignment Achieve union Promote bone growth 	<ul style="list-style-type: none"> Muscle atrophy Skin deterioration Non/mal union Joint stiffness Wound management 	<ul style="list-style-type: none"> Delta-Lite® Plus Delta-Cast® Elite Delta-Cast® Prints Delta-Cast® Soft Actimove® Supports
 <p>Supra and Sub Cast Protection</p>	<ul style="list-style-type: none"> Protect limb/soft tissue Protect cast Position limb 	<ul style="list-style-type: none"> Skin sensitivity Maceration Pressure sores Smell Cast breakdown 	<ul style="list-style-type: none"> Soffban® Natural Soffban® Synthetic Artiflex® Soft Velband™ AquaCast® Liner Delta-Dry® Protouch® Stockinette Delta Terry-Net™ S BSN Cast Shoe
 <p>Cast Removal</p>	<ul style="list-style-type: none"> Remove cast quickly and safely 	<ul style="list-style-type: none"> Saw and dust inhalation Heat build up 	<ul style="list-style-type: none"> deSoutter Cast Saws Cast Shears and Spreaders
 <p>Post Healing Rehabilitation</p>	<ul style="list-style-type: none"> Protect injury site Regain strength and mobility 	<ul style="list-style-type: none"> Protect injury site Regain strength and mobility 	<ul style="list-style-type: none"> Actimove® Genu Eco Actimove® Genu Tri-Panel Actimove® Lower Leg Walker Actimove® Umerus Comfort Actimove® Umerus Actimove® Mitella Comfort Actimove® Mitella Actimove® Mitella Arm Sling Actimove® Mitella for Kids Actimove® Sling Actimove® Sling for Kids Actimove® Vertebrace Actimove® TaloMotion Actimove® Genu Motion Actimove® ManuMotion Actimove® Manus Forte Actimove® Manus Forte Plus Actimove® EpiMotion Actimove® Cervical Actimove® Rhizo Forte Actimove® Solea Actimove® Post-Op ROM Knee Brace Actimove® Shoulder Abduction Sling

Dynacast® Prelude

Acute Synthetic Splinting System



- All in one system
- High conformability
- Easy moulding
- Smooth splint surface
- Easy to cut
- Lightweight

Dynacast® Prelude is an acute synthetic splinting system combining polyurethane coated fibreglass and polypropylene padding. Dynacast® Prelude comes pre-assembled and ready to use, giving a cost effective splinting alternative to traditional ways of providing support and immobilisation.

Single layer knit design

Dynacast® Prelude is constructed using an innovative 3D knitting process, which produces a single interlocking knit layer construction.

This construction ensures high conformability to the limb, optimum splint smoothness and increased patient comfort throughout the term of the treatment.



Indications

- Primary immobilisation of acute fractures and soft tissue injuries

Ordering information

Ref. No.	Size	Items per Unit
Roll Presentation		
71418-00	2.5cm x 4.6m	1 roll
71420-00	5cm x 4.6m	1 roll
71421-00	7.5cm x 4.6m	1 roll
71415-00	10cm x 4.6m	1 roll
71416-00	12.5cm x 4.6m	1 roll
71417-00	15cm x 4.6m	1 roll
Pre-Cut Lengths		
72869-00	5cm x 25cm	bx/5
47222-00	7.5cm x 30cm	bx/5
47225-00	10cm x 40cm	bx/5
47226-00	10cm x 76cm	bx/4
47235-00	12.5cm x 76cm	bx/4

Gypsona®

Plaster of Paris Bandage




- Fast immobilisation of fractures
- Leno cloth for good conformability and stability
- Rich creamy texture
- Low plaster loss
- Moisture resistant packaging

Gypsona® Plaster of Paris (POP) bandages are made from a leno cloth, uniformly impregnated with a high quality fast setting Plaster of Paris, which has a rich creamy texture for easy handling. The interlocked weave of the leno cloth makes the wet bandages stable, conformable and easy to mould and minimises plaster loss.

Benefits

Gypsona® retains the Plaster of Paris on the cloth so plaster loss is at a minimum whilst soaking and when being moulded. The cloth permits Gypsona® to be accurately and easily moulded to give a smooth inside surface and uniform thickness.

Gypsona® will set and harden quickly to give strong immobilisation for trauma and corrective work.



Indications

- For external splinting of fractures
- Immobilisation of diseased bones and joints
- Correction of deformities
- Emergency splinting
- Making negative and positive casts

Ordering information

Ref. No.	Size	Items per Unit
Gypsona® B.P. Plaster Bandages		
71987-51	5cm x 2.75m	12 rolls/bx
71987-54	7.5cm x 3.5m	24 rolls/bx
71987-40	10cm x 3.5m	24 rolls/bx
71987-45	15cm x 3.5m	12 rolls/bx
71987-50	20cm x 3.5m	12 rolls/bx
Gypsona® Slab Dispenser - 5 Ply		
71988-21	10cm x 20m	1 roll
71988-22	15cm x 20m	1 roll
71988-23	20cm x 20m	1 roll


Gypsona® Extra 3

Plaster of Paris Bandage



- Leno weave gauze for good conformability and stability
- Rich creamy texture
- Low plaster loss
- Moisture resistant packaging
- Rapid immobilisation of fractures

Gypsona® Extra 3 Plaster of Paris (POP) bandages are made from a leno weave gauze, uniformly impregnated with a high quality fast setting Plaster of Paris, which has a rich creamy texture for easy handling. The interlocked weave of the leno gauze makes the wet bandages stable, conformable and easy to mould as well as minimising plaster loss.



Indications

- External splinting of fractures
- Immobilisation of diseased bones and joints
- Correction of deformities
- Emergency splinting
- Making negative and positive casts

Ordering information

Ref. No.	Size	Items per Unit
71987-57	7.5cm x 3.5m	24 rolls/bx
71987-58	10cm x 3.5m	24 rolls/bx
71987-59	15cm x 3.5m	24 rolls/bx
71987-60	20cm x 3.5m	24 rolls/bx

Specialist®

Plaster of Paris Bandage



- Excellent strength
- Extra fast setting
- Low plaster loss

Specialist® Plaster of Paris uses a time-tested formula and manufacturing process to produce a consistent plaster that is uniformly coated and firmly adhered to the gauze. The result is a plaster with high dry/wet plaster retention. This means that more plaster is applied to the cast and not wasted in the wrapper or dip bucket. Specialist® casts have early green strength thus reducing the time that the cast is vulnerable to distortion

Benefits

Traditional dry formula feel and easy to use characteristics of traditional plaster.

Consistency Eliminates the need to adjust the application technique. No variation in quality from roll-to-roll.

High wet/dry plaster retention More plaster on the roll.

Early high green strength Minimises the period when casts may be vulnerable to distortion.

Strength Provides casts that are ideal for a wide range of applications.

Indications

- Fracture management
- Podiatry
- Talipes Equinovarus

Ordering information

Ref. No.	Size	Items per Unit
72258-00	7.5cm x 4.5m	12 rolls/bx
72258-01	10cm x 4.5m	12 rolls/bx
72258-02	15cm x 4.5m	12 rolls/bx
72258-03	20cm x 4.5m	12 rolls/bx

Tensoplast® Skin Traction Kit

Adhesive and Non-Adhesive Skin Traction Kits



Low Allergy, Adhesive Skin Traction Kit

- Low allergy adhesive
- Excellent conformability for easy application
- Reliable and secure
- Compression bandage for extra security

Non-Adhesive Skin Traction Kit

- High performance elastic bandage for extra security
- Good skin friction characteristics
- High conformability for easy application

Indications

- Abrasions or skin lacerations
- Impairment of circulation
- Traction weight greater than that which can be applied through the skin

Ordering information

Ref. No.	Size	Items per Unit
71544-04	Adhesive (Adult)	bx/1
71554-03	Adhesive (Child)	bx/1
71555-05	Non-Adhesive (Adult)	bx/1
71556-05	Non-Adhesive (Child)	bx/1

Delta-Lite® Plus

Fibreglass Casting Bandage



- Greater conformability
- Easy to apply
- High lamination

Delta-Lite® Plus is a fibreglass cast tape that combines a fibreglass substrate and advanced water activated polyurethane resin system which provides tack free properties. The tack free properties allow ease of application and moulding, using standard examination gloves.

Easy to Apply

The resin formulation has improved the unrolling tension of the bandage, making it easier to apply the cast, using standard examination gloves.

Greater Conformability

The heat treated substrate of Delta-Lite® Plus improves the conformability, providing greater moulding capabilities for rigid primary and secondary casting applications.

High Lamination

The resin system also ensures maximum lamination and end laydown, resulting in a smooth finished cast which is strong, long lasting and lightweight.

Indications

- Primary immobilisation of acute fractures and soft tissue injuries



Ordering information

Ref. No.	Size	Items per Unit
White		
73458-01	5cm x 3.6m	bx/10 rolls
73458-02	7.5cm x 3.6m	bx/10 rolls
73458-03	10cm x 3.6m	bx/10 rolls
73458-04	12.5cm x 3.6m	bx/10 rolls
Red		
73458-30	5cm x 3.6m	bx/10 rolls
73458-31	7.5cm x 3.6m	bx/10 rolls
73458-32	10cm x 3.6m	bx/10 rolls
Blue		
73458-20	5cm x 3.6m	bx/10 rolls
73458-21	7.5cm x 3.6m	bx/10 rolls
73458-22	10cm x 3.6m	bx/10 rolls
73458-23	12.5cm x 3.6m	bx/10 rolls
Purple		
73458-60	5cm x 3.6m	bx/10 rolls
73458-61	7.5cm x 3.6m	bx/10 rolls
73458-62	10cm x 3.6m	bx/10 rolls
Black		
73458-45	5cm x 3.6m	bx/10 rolls
73458-46	7.5cm x 3.6m	bx/10 rolls
73458-47	10cm x 3.6m	bx/10 rolls
Pink		
73458-55	5cm x 3.6m	bx/10 rolls
73458-56	7.5cm x 3.6m	bx/10 rolls
73458-57	10cm x 3.6m	bx/10 rolls
Green		
73458-25	5cm x 3.6m	bx/10 rolls
73458-26	7.5cm x 3.6m	bx/10 rolls
73458-27	10cm x 3.6m	bx/10 rolls
Mixed Pack*		
73458-50	5cm x 3.6m	bx/10 rolls
73458-51	7.5cm x 3.6m	bx/10 rolls
73458-52	10cm x 3.6m	bx/10 rolls

* Mixed pack: Contains 2 Deep Blue, 2 Light Blue, 2 Purple, 2 Deep Red, 2 Pink. Colors and availability subject to change. Please refer to the product label and/or package insert for full instructions on the safe use of these products.

Delta-Cast® Elite

Polyester Casting Bandage



- Easy to apply
- Clear radiograph viewing
- Very conformable
- Rigid support

Delta-Cast® Elite has been designed with ease of application and patient comfort in mind. Our patented fibreglass-free casting material combines strength and rigidity with soft edges for enhanced patient comfort.

Patented Fibreglass-Free Substrate

Delta-Cast® Elite has a patented polyester substrate featuring extensible yarns that easily shape to body contours, providing enhanced moulding capabilities for both primary and secondary casting applications and producing less dust during cast removal.

Patented Resin System

Ensures maximum lamination and end laydown for strong, long-lasting casts.

Indications

- Rigid primary and secondary casting applications where enhanced radiograph viewing is important



Ordering information

Ref. No.	Size	Items per Unit
White		
72279-00	5cm x 3.6m	bx/10 rolls
72279-01	7.5cm x 3.6m	bx/10 rolls
72279-02	10cm x 3.6m	bx/10 rolls
72279-03	12.5cm x 3.6m	bx/10 rolls
Red		
72275-00	5cm x 3.6m	bx/10 rolls
72275-01	7.5cm x 3.6m	bx/10 rolls
72275-02	10cm x 3.6m	bx/10 rolls
Blue		
72276-00	5cm x 3.6m	bx/10 rolls
72276-01	7.5cm x 3.6m	bx/10 rolls
72276-02	10cm x 3.6m	bx/10 rolls
Purple		
72277-00	5cm x 3.6m	bx/10 rolls
72277-01	7.5cm x 3.6m	bx/10 rolls
72277-02	10cm x 3.6m	bx/10 rolls
Black		
72278-00	5cm x 3.6m	bx/10 rolls
72278-01	7.5cm x 3.6m	bx/10 rolls
72278-02	10cm x 3.6m	bx/10 rolls
Pink		
73455-00	5cm x 3.6m	bx/10 rolls
73455-01	7.5cm x 3.6m	bx/10 rolls
73455-02	10cm x 3.6m	bx/10 rolls
Green		
73454-00	5cm x 3.6m	bx/10 rolls
73454-01	7.5cm x 3.6m	bx/10 rolls
73454-02	10cm x 3.6m	bx/10 rolls
Mixed Pack		
73457-00	5cm x 3.6m	bx/10 rolls
73457-01	7.5cm x 3.6m	bx/10 rolls
73457-02	10cm x 3.6m	bx/10 rolls

Delta-Cast® Prints

Patterned Polyester Casting Bandage



- Easy to apply
- Smooth finish, soft edges
- Clear radiograph viewing
- Available in five popular patterns

Delta-Cast® Prints has been designed with ease of application and patient comfort in mind. Our patented fibreglass-free casting material combines strength and rigidity with soft edges for enhanced patient comfort.

Delta-Cast® Prints large range of colourful patterns assist in overcoming anxiety of patients when having their cast applied.

Fibreglass-Free Substrate

Delta-Cast® Prints is comprised of polyester substrate which features extensible yarns that easily shape to body contours, providing enhanced molding capabilities for both primary and secondary casting applications and produces less dust during cast removal.

Patented Resin System

The resin formula of Delta-Cast® Prints with slip agent surfactant creates a slippery touch during working time for an optimal moulding feeling. This gives you the freedom to work with the glove of your choice, giving you a reliable precise working window.

Indications

- Rigid primary and secondary casting applications where enhanced radiograph viewing is important

Ordering information

Ref. No.		Size	Items per Unit
72273-22		Camouflage 5cm x 3.6m	bx/10 rolls
72273-06		Camouflage 7.5cm x 3.6m	bx/10 rolls
72273-20		Paw Prints 5cm x 3.6m	bx/10 rolls
72273-21		Paw Prints 7.5cm x 3.6m	bx/10 rolls
72273-07		Pastel 7.5cm x 3.6m	bx/10 rolls
72273-17		Pastel 5cm x 3.6m	bx/10 rolls
72273-02		Dinosaurs 5cm x 3.6m	bx/10 rolls
72273-05		Teddy Bears 5cm x 3.6m	bx/10 rolls

Delta-Cast® Soft

Semi-Rigid Polyester Casting Tape



- High conformability
- Controllable set time
- Superior durability
- High patient comfort
- Earlier patient mobility
- Versatile applications

Delta-Cast® Soft is fibreglass-free and was designed to provide variable support to fractures and soft tissue injuries. Its semi-rigidity allows for more flexibility and mobility, thus helping to reduce immobilisation related problems such as atrophy or joint stiffness. As a result, the overall healing process tends to be shortened.

Fibreglass-Free Substrate

Delta-Cast® Soft has multi-directional knit polyester substrate for excellent conformability, wrinkle free applications and enhanced moulding capabilities.



Indications

- Casting applications where a semi-rigid cast is required. Most commonly used for Talipes, serial casting and can be removed using scissors or by unwrapping

Ordering information

Ref. No.	Size	Items per Unit
White		
72708-02	2.5cm x 1.8m	bx/10 rolls
72708-12	5cm x 3.6m	bx/10 rolls
72708-22	7.5cm x 3.6m	bx/10 rolls
72708-32	10cm x 3.6m	bx/10 rolls
Red		
72708-08	5cm x 3.6m	bx/10 rolls
72708-18	7.5cm x 3.6m	bx/10 rolls
72708-28	10cm x 3.6m	bx/10 rolls
Blue		
72708-05	5cm x 3.6m	bx/10 rolls
72708-15	7.5cm x 3.6m	bx/10 rolls
72708-25	10cm x 3.6m	bx/10 rolls

Soffban® Natural

Natural Orthopaedic Padding



- Absorbent
- Conformable
- Easy to tear
- Soft
- Autoclavable
- Paper wrapped

Soffban® Natural has been developed to provide the conformability and easy handling normally associated with synthetic paddings, whilst providing the softness and absorbency of a natural padding. This makes it ideal for use under Gypsona® Plaster of Paris bandages.

Advantages

Absorbent – The Soffban® Natural fleece absorbs sweat and exudate, reducing the risk of skin maceration and irritation, particularly in hot climates.

Conformable – Soffban® Natural conforms well around awkward body contours, eliminating the need for pleats or tucks during application, to produce even coverage and reduce the chance of pressure points.

Easy to tear – Soffban® Natural tears and tears apart easily, aiding rapid application and reducing the risk of stressing the fracture while finishing padding application.



Soft – The 100% Rayon fleece of Soffban® Natural is soft and comfortable for the patient, providing excellent cushioning and protection.

Autoclavable – Soffban® Natural is autoclavable making it suitable for use in theatre, or for any procedures where sterility is desirable.

Paper-Wrapped – The paper wrap on Soffban® Natural does not need to be removed for autoclaving, but maintains the integrity of the roll during this process.

Indications

- Padding of orthopaedic conditions
 - POP casts
 - Synthetic casts
- Compression bandaging

Ordering information

Ref. No.	Size	Items per Unit
71472-00	5cm x 2.7m	12 rolls
71473-00	7.5cm x 2.7m	12 rolls
71469-00	10cm x 2.7m	12 rolls
71470-00	15cm x 2.7m	12 rolls
Sterile		
71477-07	7.5cm x 2.7m	ctn/112 rolls
71477-05	10cm x 2.7m	ctn/90 rolls
71477-06	15cm x 2.7m	ctn/18 rolls

Soffban® Synthetic

Water-Resistant Orthopaedic Padding



- Highly conformable for easy shaping to all body contours and joints
- Easily tearable for quick and trouble free application
- High loft for maximum cushioning and protection
- Feathers out easily, giving a smooth finish
- Low moisture absorbency minimises skin maceration

Soffban® is a comfortable, soft and smooth orthopaedic padding. It features a high loft for maximum cushioning and protection and feathers out easily, giving the padding a smooth finish.

Soffban® is highly conformable for easy shaping around awkward body contours and joints. The padding is resistant to moisture absorbency, minimising the risk of skin maceration. Soffban® is easily tearable for quick, trouble-free application and is also autoclavable, making it suitable for use in theatre or sterile procedures.



Indications

- Padding of orthopaedic conditions
 - POP casts
 - Synthetic casts
- Compression bandaging

Ordering information

Ref. No.	Size	Items per Unit
71464-01	5cm x 2.7m	12 rolls
71467-01	7.5cm x 2.7m	12 rolls
71486-05	10cm x 2.7m	12 rolls
71465-01	15cm x 2.7m	12 rolls
Sterile		
71486-07	10cm x 2.7m	ctn/90 rolls
71486-08	15cm x 2.7m	ctn/18 rolls

Artiflex® Soft

Orthopaedic Padding



- Provides excellent cushioning and soft for patient protection
- Padding allows moisture to pass through so it can evaporate
- Provides exceptional conformability and a smooth profile
- Padding can be sterilized for theatre use

Artiflex® Soft synthetic cast padding is composed of a blend of 85% polyester and 15% polyolefin providing excellent protection and comfort to the patient. Polyester retains good wicking properties, which allows it to transport moisture and rapidly dry out.



Ordering information

Ref. No.	Size	Items per Unit
73048-02	5cm x 2.7m	pk/6 rolls
73048-00	7.5cm x 2.7m	pk/6 rolls
73048-01	10cm x 2.7m	pk/6 rolls
73048-03	15cm x 2.7m	pk/6 rolls

Velband™

Absorbent Orthopaedic Padding



Velband™ is hypoallergenic bandage with a unique three-layer structure for easy application and long-term comfort. It provides excellent cushioning and moisture management under Plaster of Paris or synthetic casts, avoiding irritation and pressure sores by keeping the skin dry and comfortable.

Velband™ is easy to handle and remove, and is autoclavable for theatre use.



Indications

- Padding of orthopaedic conditions
 - POP casts
 - Synthetic casts
- Compression bandaging

Ordering information

Ref. No.	Size	Items per Unit
72315-01	7.5cm x 2.75m	12 rolls
72315-02	10cm x 2.75m	12 rolls
72315-04	15cm x 2.75m	12 rolls

AquaCast® Liner


Waterproof Cast Liner



AquaCast® is a waterproof cast liner specially designed to reduce the incidence of problems associated with synthetic casts. AquaCast® allows water to drain quickly from the ends of the cast, enabling patients to shower and bathe as normal and avoid the hassles associated with keeping traditional casts dry.



AquaCast® is highly conformable, with the adhesive coating on one side of the liner allowing for easy application and minimal slippage on the limb.



Removal of the cast is also simplified by use of the Saw Stop cut-resistant strip under the cast, which increases the safety margin during saw contact.



Indications

- Injury:
 - Non-displaced and minimally displaced fractures
 - Sprains
 - Post operative (healed incisions)
- Treatment:
 - Short-arm and below-knee casts
 - Hydrotherapy
 - Sports people and others with an active lifestyle
 - Patients facing particular care and hygiene issues, especially incontinence

Ordering information

Ref. No.	Size	Items per Unit
Waterproof Cast Liner		
73363-00	5cm x 1.7m	bx/12
73363-01	7.5cm x 1.7m	bx/12
73363-02	10cm x 1.7m	bx/12
73363-03	15cm x 1.7m	bx/12
Hipster Protective Liner		
73363-04	Size 2	pk/1
73363-05	Size 4	pk/1
73363-06	Size 6	pk/1
73363-07	Size 8	pk/1
Saw Stop Strip		
73363-08	3.1cm x 3m	1 roll
73363-09	Zip Stick	ea

Delta-Dry®

Water Resistant Cast Padding



- Conformable
- Water resistant
- Traditional cast removal with no extra components needed

The unique patented technology in Delta-Dry® has created a water resistant cast padding which allows easy application by the clinician and affords the patient all the lifestyle benefits of a cast they can get wet.

Advantages

- Allows patients to shower, swim and bathe like normal
- No special tools are needed for application or removal
- Conformable enabling patient to having increased comfort and protection
- Easy to apply
- Better long term water resistance
- Offers more under cast protection
- Eliminates the need for stockinette



Indications

- Non-displaced and minimally displaced fractures
- Sprains
- Post operative (healed incisions)

Ordering information

Ref. No.	Size	Items per Unit
Bandage		
73443-03	5cm x 2.4m	bx/12
73443-04	7.5cm x 2.4m	bx/12
73443-02	10cm x 2.4m	bx/12

Softliner		
72779-01	2.5cm x 10m	bx/1
72779-02	5cm x 10m	bx/1
72779-03	7.5cm x 10m	bx/1
72779-04	10cm x 10m	bx/1
72779-05	12.5cm x 10m	bx/1

Stockinette		
74564-00	2.5cm x 10m	bx/1
74564-01	5cm x 10m	bx/1
74564-02	7.5cm x 10m	bx/1
74564-03	10cm x 10m	bx/1
74564-04	12.5cm x 10m	bx/1

Protouch® Stockinette

Synthetic and Cotton Stockinette



- Available in synthetic or cotton
- Conformable and easy to use
- Available in a wide range of sizes

Protouch® synthetic and cotton stockinettes provide patients with optimum comfort at the critical juncture where the skin meets the cast.

Protouch® Synthetic Stockinette

- Made of synthetic fibres to retract moisture away from the skin for rapid evaporation
- Soft and conformable providing good cushioning and comfort to the patient
- Stretches to three times its normal width to conform to an arm, leg or body
- Available in a variety of sizes for all casting applications

Protouch® Cotton Stockinette

- Ideal material for an initial undercast layer
- Made of a high-grade knitted cotton
- Soft and conformable providing good cushioning and comfort to the patient
- Stretches to three times its normal width to conform to an arm, leg or body
- Available in a variety of sizes for all casting applications

Ordering information

Ref. No.	Size	Items per Unit
Synthetic Stockinette		
72026-00	2.5cm x 22.75m	bx/1
72026-01	5cm x 22.75m	bx/1
72026-02	7.5cm x 22.75m	bx/1
72026-03	10cm x 22.75m	bx/1
72026-05	20cm x 22.75m	bx/1
Cotton Stockinette		
72022-00	2.5cm x 22.75m	bx/1
72022-01	5cm x 22.75m	bx/1
72022-02	7.5cm x 22.75m	bx/1
72022-03	10cm x 22.75m	bx/1

Delta Terry-Net™ S

Terry Cloth Stockinette



- Conformable and easy to apply
- Cost-effective roll presentation
- Soft for patient comfort and protection
- Lightweight, flexible and durable

Delta Terry-Net™ S is for use where extra patient comfort is required. Supplied in an easy-to-use dispensing box for improved storage and handling, Delta Terry-Net™ S is latex free.

Terry Knit Design

The terry knit design of Delta Terry-Net™ S allows clinicians greater flexibility to accommodate varying limb sizes as well as increased treatment flexibility in the management of specialised casting indications, such as total contact casting. The terry knit design also allows for enhanced protection to fragile skin when managing skin tears.

Terry Knit Fibres

The soft moisture wicking fibres of Delta Terry-Net™ S provide improved patient comfort and reduced maceration whilst providing increased protection both in-situ and upon removal.

Ordering information

Ref. No.	Size	Items per Unit
72306-00	5cm x 13.5m	bx/1
72306-01	7.5cm x 13.5m	bx/1
72306-02	10cm x 13.5m	bx/1

deSoutter

Cast Removal System – CC4



The powerful CC4 CleanCast™ saw and dust extraction system comprises of the CC4 mains voltage oscillating saw, extraction hood and the CCS portable extractor unit with highly efficient dust filtration capability.

Powerful 110 Watt Motor

Maintains rapid and efficient cutting on the thickest of body casts

Control and Maneuverability

Knurled grip close to saw blade provides operator with greater flexibility of movement.

Visibility

Clear polycarbonate extraction hood offers excellent visibility when cutting. The hood is fully rotatable to suit the most convenient cutting position.

Low Noise

The CCS Extractor Unit is exceptionally quiet at only 58dB(A) and has two speed settings.

Four Levels of Filtration

A filtration level of 99.997% on all dust particles greater than 0.5 microns is achieved with the extractor.

CC4 Saw Specifications

Oscillating speed: 14000 osc/min.
 Voltage: 230V
 Weight: 1.3 kg
 Sound level: 70 dB(A)
 Length: 333 mm



Ordering information

Ref. No.	deSoutter Code	Description	Items per Unit
CC4 Saw			
72030-06	1281444	Cast Saw deSoutter Model CC4	1
72030-05	CC4A-230AUS	Cast Saw deSoutter Model CC4 (Including Extractor Unit)	1
72035-25	9490AUS	Dust Extractor Unit Model 230v	1
72035-34	78472	Dust Extractor Hose Kit for CC4	1
72035-01	602763	Dust Bags for CCS Extractor	1
72033-08	6250	Cast Saw Blade to suit deSoutter 64mm	1

deSoutter

Cast Removal System – CC5



The CleanCast™ low voltage saw and integrated extractor system combines state-of-the-art ergonomics with the performance and durability you’ve come to expect with all deSoutter Medical products.

Low Voltage for Increased Safety

Saw operates at only 25v d.c.

Integral Dust Extraction Path

Vacuum is ducted internally through the handpiece to provide better visibility and handling.

Two Speed Control on Both Saw and Extractor

Quiet, low speed setting is ideal for use with children, whilst high speed setting allows thick body casts to be cut quickly.

Four Levels of Filtration

A filtration level of 99.997% on all dust particles greater than 0.5 microns is achieved with the extractor.

Low Noise

Specially designed rubber motor mounts insulate the internal mechanism from the outer case. This results in an extremely low sound.

CC5 Saw Specifications

High Oscillating speed: 16000 osc/min.
 Low Oscillating speed: 11800 osc/min.
 Voltage: 230V
 Weight: 0.85 kg
 Sound Level: 68 dB(A)
 Length: 245 mm



Ordering information

Ref. No.	deSoutter Code	Description	Items per Unit
CC5 Saw			
72036-10	CC5A-230AUS	Cast Saw deSoutter Model CC5 (Including Extractor Unit)	1
72035-25	9490AUS	Dust Extractor Unit Model 230v	1
72034-00	6400	Dust Extractor Hose Kit for CC5	1
72035-01	602763	Dust Bags for CCS Extractor	1
72033-08	6250	Cast Saw Blade to suit deSoutter 64mm	1

deSoutter

Battery Powered Saw – CSB-100



The deSoutter Medical CSB-100 saw permits the user total freedom of movement which is particularly useful for emergency departments, ward work, veterinary clinics and numerous other applications where the cordless concept is the only practical solution.

CSB-100 Saw Specifications

Oscillating speed: 16000 osc/min.
Voltage: 10.8 Vdc
Weight: 1.05 kg

Battery Charger Specifications

Supply voltage: 100-240 Vac
Supplier frequency: 50-60 Hz


Supplied with interchangeable mains supply adaptors for use in most countries.

Ordering information


Ref. No.	deSoutter Code
CSB-100 Cast Saw	
1286774	73065-00
Includes:	CSB-100 Saw
	CB-100 Battery x 2
	BC-100 charger
	Saw Blade
	Hexagon Wrench
	Carry Case

deSoutter

Battery/Mains Powered Saw – CSP-201



25mm
Ref No. 15310
73356-01



29mm
Ref No. 15320
73356-02



Specifically developed for paediatric use, the Red Dot award winning CSP System is exceptionally quiet and compact minimising patient anxiety. The CSP System offers both battery and mains power options, and lithium ion battery technology provides an unrivalled power to weight ratio.

CSP-1201 Saw Specifications

Oscillating speed: 14800-17500 osc/min.
Weight: 319 grams
Sound level: 55 dB(A)

Battery Specifications

Capacity: 30.8 Wh
Weight: 437 grams
Charging time from flat: < 3 hours

Supplied with interchangeable mains supply adaptors for use in most countries.

- Soft start feature with noise control helps minimise patient anxiety
- Continuous 2 hour – cutting performance
- Built-in service reminder – 300 charging cycles
- Display allows the user to select the optimum cutting speed and view the remaining battery capacity
- Micro processor controlled voltage booster ensures constant cutting performance - zero deterioration

Ordering information

Ref. No.	deSoutter Code
CSP-201 Pro Cast Saw	
1286994	73356-00
Includes:	CSP-201 Saw and Cable
	15250 BP-201 Battery Pack Pro
	15270 IC-201 Inter Connector
	15280 PS-201 Power Supply
	15290 Carry Case
	15300 Wrench
	73356-01 15310 Saw Blade
	73356-02 15320 Saw Blade

Cast Room Accessories

Cast Saw Blades

- Available in PTFE (Teflon), extra durable titanium nitride and stainless steel coatings
- Range of sizes and profiles to cover all requirements
- Suits deSoutter Cast Saws



Cast Spreaders

- High quality surgical stainless steel construction
- Spring loaded handles for easy use

Cast Removal Shears

- Ideal for cutting casting materials, paddings, felt, etc
- High quality surgical stainless steel construction



Cast Scissors

- Lister 200mm bandage scissors
- High quality surgical stainless steel construction

Ordering information

Ref. No.			Items per Unit
Cast Saw Blades			
6230	Circular PTFE Coated	44mm	1
72033-25 6240	Circular PTFE Coated	50mm	1
72033-08 6250	Circular PTFE Coated	64mm	1
72033-03 6260	Segment PTFE Coated	80mm	1
72033-10 67400065	Circular Stainless Steel	64mm	1
Cast Spreaders			
72046-11 5120		Small	1
73066-00 5130		Large	1
Cast Removal Shears			
73066-01 7680			1
Cast Scissors			
72937-00		20cm Lister	1
72294-02		20cm Clean Cut	1

BSN Cast Shoe

Open Toe Cast Shoe



- Durable
- Precise fit
- Rocker bottom

Durable Cast Shoes

The BSN Cast Shoe is constructed with quality materials using the highest manufacturing standards. Constructed with a heavy canvas upper, it is a premium cast shoe at an affordable price. Three hook and loop enclosures ensure a precise fit while the rocker bottom helps promote normal gait pattern.

Advantages

- High level of durability
- Precise fit
- Maximum shock absorption
- Comfortable and easy to apply

Ordering information

Ref. No.		Size	Items per Unit
71427- 23	Toddler	16.5cm x 9cm	1 shoe
71427- 24	Extra Small	19cm x 10cm	1 shoe
71427- 25	Small	22cm x 10.5cm	1 shoe
71427- 26	Medium	26.5cm x 11.5cm	1 shoe
71427- 27	Large	30cm x 13cm	1 shoe
71427- 28	Extra Large	32cm x 13.5cm	1 shoe

Cast Protection

Actimove® Post-Op ROM Knee Brace

Knee Brace



- Easy to apply thanks to the hook and loop buckle closure system
- Provides comfort to the patient, as it is light weight and made from breathable materials
- Latex-free formula
- The brace can be universally worn, left or right

Actimove® Post-Op ROM Knee Brace is designed to provide secure, comfortable knee immobilisation or range of motion limitation for a range of patients through its universal fit and adjustability to various leg lengths.

Lightweight, it optimally immobilises the patient during their treatment, suitable for use post-injury, post-operatively, or throughout rehabilitation.

Advantages

- For a faster application, the hinge offers quick lock keys to keep the knee at 0°, 15° and 30° degrees
- To limit knee movement more precisely, the range of motion can be set in flexion from -10° to 120° and extension from -10° to 90° in increments of 10°
- The hinge bar locking system provides security to users
- The telescopic hinge bars of the brace can be extended from 44.5cm - 66 cm to fit most patients

Indications

For post-operative use after knee procedures like:

- ACL, PCL, MCL, LCL injuries
- Tibial plateau fractures
- Osteochondral repairs
- Meniscal repairs
- Patella tendon repairs
- Condylar fractures
- Sprains or strains of the knee
- HTOs (High Tibial Osteotomy)

Ordering information

Ref. No.	Items per Unit
76278-55	1 knee brace

Actimove® Genu

ECO Knee Immobiliser



Orthopaedic Soft Goods

- Comfortable and skin-friendly
- Low profile and effective patient fit
- Washable

Actimove® Genu ECO Knee Immobiliser supports the knee and helps to reduce pain after trauma or surgery.

Advantages

- Pre-shaped foam for knee and calf
- A detachable hook and loop sensitive closure system with D-rings permits individual patient adjustments
- Removable and adjustable aluminium stays allow good anatomical fit
- Terry cloth interior padding for patient comfort
- Splint is washable for additional patient hygiene


Extension Tabs

- Foam tabs, with hook, attach to ECO Knee Immobiliser to increase the width of the device
- This system fits wider leg circumferences, reducing inventory of sizes

Indications

- Following surgery, soft tissue injuries and minor fractures in the knee region

Ordering information Sizing Chart



Ref. No.	Size	Measurements
72819-10	B1-S	< 180cm 60cm < 52cm
72819-11	B2-L	< 180cm 60cm 45-65cm
72819-12	C1-S	> 180cm 70cm < 52cm
72819-13	C2-L	> 180cm 70cm 45-65cm
72819-14		50cm < 64cm
72819-15		50cm > 64cm
72819-16		60cm < 64cm
72819-17		60cm > 64cm

Eco Knee Extension Tabs

72819-74	Small	pk/10
72819-75	Medium	pk/10
72819-76	Large	pk/10

*Extension Tabs exist in three sizes adapted to the three lengths of the Actimove ECO Knee Immobiliser.

Actimove® Genu

Tri-Panel Knee Immobiliser



- Two sizes fit all
- Effective patient fit
- Washable

Actimove® Genu Tri-Panel Knee Immobiliser supports the knee and helps to reduce pain after trauma or surgery.


Special Design Features

- Panels allow splint adjustments for a broad spectrum of patient leg circumferences
- Hook and loop sensitive closure system permits easy application and removal
- Removable and adjustable aluminium stays allow good anatomical fit
- Terry cloth interior padding for patient comfort
- Splint is washable for additional patient hygiene

Indications

- Following surgery, soft tissue injuries and minor fractures in the knee region

Ordering information Sizing Chart



Ref. No.	Size	Measurements
72819-14	A1-S 	< 160cm 50cm < 62.5cm
72819-15	A2-L 	< 160cm 50cm > 62.5cm
72819-16	B1-S 	< 160cm 60cm < 62.5cm
72819-17	B2-L 	< 160cm 60cm > 62.5cm

Actimove®

Lower Leg Walker



- Quick and easy application
- Superior comfort

A high-quality rigid outer skeleton ensures high protection of the foot, ankle and lower leg. Circumferential strapping, a shock absorbing sole and breathable materials provide great comfort for the patient. The rocker bottom of the Actimove® Walker facilitates a natural physiological gait when walking.

Actimove® Air Walker

Actimove® Air Walker has a liner which can be inflated for oedema control. In just a few pumps, the inflated liner can reduce the void around the foot and ankle. This improves intimacy of fit, reduces unwanted movement, and allows for optimal patient comfort and wearing tolerance.

Comfort

Anatomically contoured struts prevent ankle motion without causing friction against the sides of the calf. A full-length insole of pressure mouldable comfort foam conforms to the foot.

Indications

- Available in 5 sizes (XS, S, M, L, XL)
- Suited for a broad range of post operative and post trauma applications

Ordering information

Ref.No.	Size	UK	US	EU
Standard Walker				
High	Low	Shoe Size		
76270-35	76270-30	XS	< 3.5 < 3	< 4.5 < 6 < 36.5 < 36
76270-36	76270-31	S	4-6.5 3.5-5.5	4.5-7 6-8 36-39.5 36.5-38.5
76270-37	76270-32	M	7-9.5 6-8.5	7.5-10 8.5-11 40-44 39-43
76270-38	76270-33	L	10-12 9-11	10.5-12.5 11.5-13.5 44.5-46.5 44-46
76270-39	76270-34	XL	> 12 > 11	> 12.5 > 13.5 > 46.5 > 46
Air Walker				
High	Low	Shoe Size		
76272-35	76272-30	XS	< 3.5 < 3	< 4.5 < 6 < 36.5 < 36
76272-36	76272-31	S	4-6.5 3.5-5.5	4.5-7 6-8 36-39.5 36.5-38.5
76272-37	76272-32	M	7-9.5 6-8.5	7.5-10 8.5-11 40-44 39-43
76272-38	76272-33	L	10-12 9-11	10.5-12.5 11.5-13.5 44.5-46.5 44-46
76272-39	76272-34	XL	> 12 > 11	> 12.5 > 13.5 > 46.5 > 46

Actimove® Solea

Post-Op Shoe



- Wide construction allows application of bulky bandage or wound dressing
- Rigid sole reduces foot and toe flexion
- Soft and breathable liner
- Extra padding in the heel area for advanced comfort
- The hook and loop closure system facilitates easy application
- Anti-slip shoe sole provides a safe surface for patient ambulation
- Universal fit for left and right foot

Indications

- To help protect the foot following surgery of the foot and/or/toes

Ordering information

Ref. No.	Size	EU	UK	Items per Unit
		Shoe Size		
72819-51	1	36	3	pk/1
72819-52	2	37-38	5	pk/1
72819-53	3	39-40	6-7	pk/1
72819-54	4	41-42	8	pk/1
72819-55	5	43-44	9	pk/1
72819-56	6	45	10	pk/1

Actimove® Shoulder Abduction Sling

Shoulder Immobiliser



- Easy to apply and to remove, even with one hand, thanks to the stable snap-lock buckle connecting the strap with the pocket, and the quick release buckle attached to the waist strap
- Designed to aid pain relief through shoulder immobilisation at the upper-arm abduction at around 15°
- Secure and tailored fit for most patients through the extra-long waist strap of which can be adjusted in length
- Increased comfort through the adjustable padded neck strap reducing the pressure on the neck
- For increased patient comfort during prolonged wear as the perforated material allows for breathability
- Brings comfort during wear due to the soft outer and inner materials, edges and seams
- Helps to maintain muscle strength through hand and arm muscle training exercises using the attached squeeze ball
- Clothes are protected through the soft, supple and colour fast material
- Dirt and signs of wear are less visible due to the black inner lining of the arm pocket

Actimove® Shoulder Abduction Sling is designed to aid pain relief by immobilizing the shoulder joint. It limits the range of motion, providing rest for the shoulder in a flexed elbow position and upper-arm abduction at around 15°.

Patients will appreciate the easy one-handed application and the comfort of the soft, breathable and washable material.

Indications

- Minor instability of the shoulder (e.g., dorsal instability)
- Postoperative use (e.g., after rotator cuff rupture)
- Joint dislocations (e.g., shoulder luxation)



Easy one-handed application thanks to the stable snap-lock buckle.

**LATEX
FREE
FORMULA**

Ordering information

Sizing Chart



Ref. No.	Size	Measurements	Items per Unit
73446-04	S	< 37cm	1 immobiliser
73446-05	M	37–40cm	1 immobiliser
73446-06	L	40–43cm	1 immobiliser
73446-07	XL	> 43cm	1 immobiliser

Actimove® Umerus Comfort

Shoulder Immobiliser



- Helps to prevent patients from experiencing neck pressure and the hassle of a moving neck pad thanks to the broad 5.5 cm wide sling
- Improves shoulder stabilisation through the waist strap which can be adjusted to suit most patients
- Allows for better hand positioning through the thumb loop which limits the rotation of the forearm while still allowing full finger movement
- Simple to apply and adjust for an individual fit with easy threading through the freely rotating O-rings
- Patients with a compromised arm can easily take the immobiliser off and back on, as a one-handed self-application is helped by the hook and loop Y-Tabs and open-and-close pocket design
- Cares for clothes through the gentle pocket material
- Machine washable at 40°C
- Radiolucent
- Metal-free

Actimove® Umerus Comfort is the superior comfort shoulder and upper arm immobiliser for limited abduction, extension, internal and external rotation of the humerus. Patients will appreciate the easy application and breathable material while they are immobilised.

Indications

- Postoperative support
- Fracture management (e.g. clavicle, scapular, humeral, and greater tuberosity fractures)
- Joint dislocations (e.g. acromioclavicular & sternoclavicular dislocations, shoulder luxation)
- Inflammatory diseases (e.g. calcific tendinitis of the rotator cuff)
- Minor instability of the shoulder
- Support for hemiplegic patients

**LATEX
FREE
FORMULA**

Ordering information

Sizing Chart



Ref. No.	Size	Measurements	Items per Unit
72819-94	XS	< 33cm	1 immobiliser
72819-39	S	< 37cm	1 immobiliser
72819-40	M	37-40cm	1 immobiliser
72819-41	L	40-43cm	1 immobiliser
72819-42	XL	> 43cm	1 immobiliser

Actimove® Umerus

Shoulder Immobiliser



- Improves shoulder stabilisation through the waist strap which can be adjusted to suit most patients
- Allows for better hand positioning thanks to the thumb loop which limits the rotation of the forearm while still allowing full finger movement
- Simple to apply and adjust for an individual fit with easy threading through the freely rotating O-rings
- Machine washable at 40°C
- Radiolucent
- Metal-free

Actimove® Umerus is the cost-effective shoulder and arm immobiliser for limited abduction, extension, internal and external rotation of the humerus.

It helps to keep your inventory low, as it fits on both left and right arm, available in four sizes for your patients who will appreciate its easy application.



Indications

- Post operative support
- Fracture management (e.g. clavicle, scapular, humeral, and greater tuberosity fractures)
- Joint dislocations (e.g. acromioclavicular & sternoclavicular dislocations, shoulder luxation)
- Inflammatory diseases (e.g. calcific tendinitis of the rotator cuff)
- Minor instability of the shoulder
- Support for hemiplegic patients

**LATEX
FREE
FORMULA**

Ordering information

Sizing Chart



Ref. No.	Size	Measurements	Items per Unit
72819-65	S	< 37cm	1 immobiliser
72819-66	M	37-40cm	1 immobiliser
72819-67	L	40-43cm	1 immobiliser
72819-68	XL	> 43cm	1 immobiliser

Actimove® Mitella Comfort

Arm Sling



- Helps to prevent patients from experiencing neck pressure and the hassle of a moving neck pad thanks to the broad 5.5cm wide sling
- Allows for better hand positioning through the thumb loop which limits the rotation of the forearm while still allowing full finger movement
- Simple to apply and adjust for an individual fit with easy threading through the freely rotating O-rings
- Patients with a compromised arm can easily take the support off and back on, as a one-handed self-application is helped by the hook and loop Y-Tabs and open-and-close pocket design
- Cares for clothes through the gentle pocket material, while also protecting casts or bandaging from moisture and dirt
- Machine washable at 40°C
- Radiolucent
- Metal-free

Actimove® Mitella Comfort is the superior comfort arm support. Patients will appreciate the easy application and breathable material.

Indications

- Postoperative support
- Fracture management (e.g. clavicle, scapular, humeral, and greater tuberosity fractures)
- Joint dislocations (e.g. acromioclavicular & sternoclavicular dislocations, shoulder luxation)
- Inflammatory diseases (e.g. calcific tendinitis of the rotator cuff)
- Minor instability of the shoulder
- Support for hemiplegic patients
- Support while wearing a cast or splint
- Temporary arm stabilisation in emergency/paramedic care

**LATEX
FREE
FORMULA**

Ordering information

Sizing Chart



Ref. No.	Size	Measurements	Items per Unit
72819-61	S	< 37cm	1 support
72819-62	M	37-40cm	1 support
72819-63	L	40-43cm	1 support
72819-64	XL	> 43cm	1 support

Actimove® Mitella

Arm Sling



- Allows for better hand positioning thanks to the thumb loop which limits the rotation of the forearm while still allowing full finger movement
- Simple to apply and adjust for an individual fit with easy threading through the freely rotating O-rings
- Cares for clothes through the gentle pocket material, while also protecting casts or bandaging from moisture and dirt
- Machine washable at 40°C
- Radiolucent
- Metal-free

Actimove® Mitella is the cost-effective arm support. It helps to keep your inventory low, as it fits on both left and right arm, available in four sizes for your patients who will appreciate its easy application.

Indications

- Postoperative support
- Fracture management (e.g. clavicle, scapular, humeral, and greater tuberosity fractures)
- Joint dislocations (e.g. acromioclavicular and sternoclavicular dislocations, shoulder luxation)
- Inflammatory diseases (e.g. calcific tendinitis of the rotator cuff)
- Minor instability of the shoulder
- Support for hemiplegic patients
- Support while wearing a cast or splint
- Temporary arm stabilisation in emergency/paramedic care

**LATEX
FREE
FORMULA**

Ordering information

Sizing Chart



Ref. No.	Size	Measurements	Items per Unit
72819-57	S	< 37cm	1 immobiliser
72819-58	M	37-40cm	1 immobiliser
72819-59	L	40-43cm	1 immobiliser
72819-60	XL	> 43cm	1 immobiliser

Actimove® Mitella for Kids

Arm Sling



- Distracts girls and boys from an often traumatic situation through the unique, playful pattern
- Helps to drive compliance in children by providing comfort while wearing
- Simple to apply and adjust for an individual fit with easy threading through the freely rotating O-rings
- Cares for clothes through the gentle pocket material, while also protecting casts or bandaging from moisture and dirt
- Machine washable at 40°C
- Radiolucent
- Metal-free

Actimove® Mitella for Kids is the child-friendly arm support with a fun and colourful pocket to make kids, parents and clinicians smile. It provides stable support and fits on both left and right arm with universal paediatric sizing, fitting most smaller patients.

Indications

- Support while wearing a cast or splint
- Fracture management (e.g. clavicle, scapular, humeral, and greater tuberosity fractures)
- Temporary arm stabilisation in emergency/paramedic care
- Joint dislocations (e.g. acromioclavicular and sternoclavicular dislocations, shoulder luxation)
- Postoperative support
- Minor instability of the shoulder
- Support for hemiplegic patients

**LATEX
FREE
FORMULA**

Ordering information

Ref. No.	Size	Items per Unit
74561-13	One size - child	1 support

Actimove® Sling

Collar and Cuff Support



- Comfortable and secure support
- Strong, secure “Y” closure tabs
- Convenient dispenser box presentation

The padded strap and “Y” closure tabs of the Actimove® Sling provide effective vertical support to the collarbone and upper limb area.

Advantages:

- Extra patient comfort with padded strap
- Y-Tab hook system facilitates secure clinical positioning
- Convenient dispenser box presentation
- 28 Y-Tabs are included in each box. Additional Y-Tabs can be purchased, if required



Indications

Following trauma injuries, collarbone injuries and soft tissue injuries

Ordering information

Ref. No.	Size	Items per Unit
72859-18	2 rolls x 12m	2 rolls/bx
Y-Tabs		
73466-00	14 tabs/pkt	pk/1

Actimove® Sling for Kids

Collar and Cuff Support



- Distracts girls and boys from an often traumatic situation through the unique, playful pattern
- Easily remove the sling from the poly bag by cutting an opening on top for access to the middle of the roll
- Drives patient compliance through increased comfort, as rolling of the sling is prevented by the three lengthwise seams
- When removing the sling, this can be done without opening the Y-Tabs (includes 14 Y-Tabs)
- Machine washable at 40°C
- Radiolucent
- Metal-free

Actimove® Sling for Kids is the child-friendly arm sling in a fun and colourful pattern to make kids, parents and clinicians smile.

The one-roll polybag pack uses little storage space and allows for easy access to the sling through its top opening system. Individually adjustable in length, provide arm slings for a stable and supported forearm through various application styles in one product.

Indications

- Fracture management (e.g. clavicle, scapular, humeral, and greater tuberosity fractures)
- Support while wearing a cast or splint
- Joint dislocations (e.g. acromioclavicular & sternoclavicular dislocations, shoulder luxation)
- Temporary arm stabilisation in emergency/paramedic care

**LATEX
FREE
FORMULA**

Ordering information

Ref. No.	Size	Items per Unit
73117-00	3.6cm x 10.8m	1 roll in a dispenser bag

Actimove® TaloMotion

Functional Ankle Support



- Breathable knitting structure for better patient comfort
- 3D anatomical fit for improved wearing comfort
- Special knitting structure in the sole area for higher stability against abrasion
- Seamless integration of pelottes, lined with soft terry cloth
- Soft latex-free materials for enhanced comfort
- Attractive design and colour combination



Advanced knitting structure provides gradient compression and increased stability



Two anatomically shaped pelottes provide targeted pressure to reduce swelling and pain



Pressure-free cuffs without elastic band prevent constriction



Indications

- Postoperative and post traumatic irritation
- Joint effusion and swelling
- Ankle arthrosis and arthritis
- Tendomyopathy
- Slight ligament instabilities

**LATEX
FREE
FORMULA**

Ordering information

Sizing Chart



Ref. No.	Size	Measurements	Items per Unit
Left Foot			
73487-42	XS	17–19cm	pk/1
73487-43	S	19–21cm	pk/1
73487-44	M	21–23cm	pk/1
73487-45	L	23–25cm	pk/1
73487-46	XL	25–27cm	pk/1
Right Foot			
73487-12	XS	17–19cm	pk/1
73487-13	S	19–21cm	pk/1
73487-14	M	21–23cm	pk/1
73487-15	L	23–25cm	pk/1
73487-16	XL	25–27cm	pk/1

Actimove® GenuMotion

Functional Knee Support



- Innovative, award-winning Helix Design
- 3D anatomic fit for better fit and wearing comfort
- Anatomically shaped pelotte stabilises the patella and improves proprioception
- Spiral stays support lateral stabilisation
- Seamless integration of stays and pelotte, lined with soft terry cloth
- Breathable knitting structure for better patient comfort
- Soft, latex-free materials for enhanced comfort


Indications

- Arthrosis
- Arthritis
- Minor instabilities of the knee joint
- Minor meniscus damages
- Weak ligaments

**LATEX
FREE
FORMULA**



Advanced knitting structure provides gradient compression supporting pain relief



Specially designed Motion Comfort Zone in the popliteal area, pleasantly soft and almost pressure free



Pressure-free cuffs for a non-slip fit with no constriction

Ordering information Sizing Chart



Ref. No.	Size	Measurements	Items per Unit
73468-11	XS	38-41cm 28-31cm	pk/1
73468-12	S	41-44cm 31-34cm	pk/1
73468-13	M	44-47cm 34-37cm	pk/1
73468-14	L	47-50cm 37-40cm	pk/1
73468-15	XL	50-53cm 40-43cm	pk/1
73468-16	XXL	53-56cm 43-46cm	pk/1

Actimove® ManuMotion

Functional Wrist Support



Orthopaedic Soft Goods

- Advanced knitted structure for targeted compression
- Mouldable stay stabilises the wrist in a functional position and supports pain relief
- Dorsal pressure pad supports the natural lymph flow and helps reduce swelling
- Soft motion comfort zone in the thumb area for high wearing comfort
- Inelastic wrist band for additional wrist stabilisation

Indications

- Arthritis, rheumatoid arthritis
- Tenosynovitis
- Tendomyopathy, Tendopathy
- Post operative
- Irritations of soft tissue
- Distortion, contusion

**LATEX
FREE
FORMULA**

Ordering information

Sizing Chart



Measure wrist circumference at wrist bone

Ref. No.	Size	Measurements	Items per Unit
73497-12	Right XS	14-15cm	pk/1
73497-13	Right S	15-16cm	pk/1
73497-14	Right M	16-17cm	pk/1
73497-15	Right L	17-18cm	pk/1
73497-16	Right XL	18-19cm	pk/1
73497-42	Left XS	14-15cm	pk/1
73497-43	Left S	15-16cm	pk/1
73497-44	Left M	16-17cm	pk/1
73497-45	Left L	17-18cm	pk/1
73497-46	Left XL	18-19cm	pk/1

Actimove®

Rhizo Forte



- Firm cervical support
- Easy to apply
- Convenient colour coded identification system

Actimove® Rhizo Forte brace is designed to restrict thumb movement whilst maintaining the daily functionality of the hand. The Red Dot® award winning design ensures the immobilisation of the CMC and MCP joints without restricting the movement of the wrist and the other four fingers.

Special features:

- Low profile design with smooth edges for high wearing comfort
- Composition and design that is easy to clean
- Closures which allow the brace to be easily secured with one hand
- Effectively restricts thumb movement, while leaving the four fingers and the wrist to move freely
- Mouldable aluminium core is easily adjustable for a perfect fit



reddot award 2014
winner

Indications

- Painful conditions and irritations of the thumb joints (e.g. Osteoarthritis)
- Rheumatoid Arthritis of the thumb
- Ulnar Collateral Ligament injuries of the MCP1 Joint (e.g. Skier's thumb)
- Soft tissue injuries (e.g. sprains)
- Ligament instability
- Post operative immobilisation

Ordering information

Sizing Chart



Ref. No.	Size	Measurements	Items per Unit
76238-00	S Right	4.5cm–5.4cm	bx/1
76238-01	S Left	4.5cm–5.4cm	bx/1
76238-02	M Right	5.5cm–6.3cm	bx/1
76238-03	M Left	5.5cm–6.3cm	bx/1
76238-04	L Right	6.4cm–7.3cm	bx/1
76238-05	L Left	6.4cm–7.3cm	bx/1

Actimove® Manus Forte

Wrist Immobiliser



Orthopaedic Soft Goods

- Secure wrist stabilisation supports pain relief
- Anatomically shaped aluminium stay, mouldable for individual fit
- Wide thumb opening with soft edges for high wearing comfort
- Smart slip-on construction, easy to apply with one hand
- Security of hook & loop closures prevent straps from slipping out of their buckles
- Soft, breathable materials to support the natural skin climate

Indications

- Pain and irritation of the wrist joint (i.e. carpal tunnel syndrome)
- Post operative and post traumatic irritation
- Tendinitis, Tenosynovitis
- Osteoarthritis of the wrist

**LATEX
FREE
FORMULA**

Ordering information Sizing Chart



Measure wrist circumference at wrist bone

Ref. No.	Size	Measurements	Items per Unit
73482-01	S Right	≤ 14cm	bx/1
73482-02	S Left	≤ 14cm	bx/1
73482-03	S/M Right	14-19cm	bx/1
73482-04	S/M Left	14-19cm	bx/1
73482-05	L/XL Right	19-23cm	bx/1
73482-06	L/XL Left	19-23cm	bx/1

Actimove® Manus Forte Plus

Wrist & Thumb Immobiliser



- Secure wrist and thumb stabilisation supports pain relief
- Two anatomically shaped aluminium stays, mouldable for individual fit
- Flat thumb seam for high wearing comfort
- Inelastic wrist strap runs counter-clockwise to prevent rotation of the brace
- Special security hook & loop closures prevent straps from slipping out of their buckles

Indications

- Irritation of the thumb saddle joint and/or proximal MCP joint of the thumb
- Acute and chronic pain and irritation in wrist and thumb joint
- Arthrosis
- Tenosynovitis
- Post operative and post traumatic irritation

**LATEX
FREE
FORMULA**



Smart slip-on construction, easy to apply with one hand



Soft, breathable materials to support the natural skin climate



Adjustable diagonal thumb strap for effective stabilisation

Ordering information Sizing Chart



Measure wrist circumference at wrist bone

Ref. No.	Size	Measurements	Items per Unit
73496-03	S Right	≤ 14cm	bx/1
73496-33	S Left	≤ 14cm	bx/1
73496-04	S/M Right	14-19cm	bx/1
73496-34	S/M Left	14-19cm	bx/1
73496-06	L/XL Right	19-23cm	bx/1
73496-36	L/XL Left	19-23cm	bx/1

Actimove® EpiMotion

Functional Elbow Support



Orthopaedic Soft Goods

- Supports pain relief, helps reduce swelling and promotes the healing process
- Two anatomically shaped pelottes provide targeted pressure to the epicondylar tendons
- Seamless terry cloth liner for enhanced comfort
- Epicondyle cut-outs and helix design facilitate correct positioning
- High-quality latex-free materials for breathability and soft touch
- Attractive design and color combination improve patient compliance

Indications

- Epicondylopathies, e.g. tennis or golfer's elbow
- Degenerative changes in the elbow
- Joint inflammation
- Post operative or post traumatic irritation
- Tendomyopathies

**LATEX
FREE
FORMULA**



Advanced knitted structure providing gradient compression



Soft Motion Comfort Zone in the elbow crease for a high degree of wearing comfort



Pressure-free cuffs without elastic band prevent constriction

Ordering information Sizing Chart



Ref. No.	Size	Measurements	Items per Unit
73477-12	XS	19-21cm	pk/1
73477-13	S	21-23cm	pk/1
73477-14	M	23-25cm	pk/1
73477-15	L	25-27cm	pk/1
73477-16	XL	27-29cm	pk/1
73477-17	XXL	29-31cm	pk/1

Actimove® Cervical

Cervical Collar



- Pre-shaped foam for anatomical fit
- Soft cover material for patient comfort
- Covers a wide range of neck sizes
- Available in medium density for medium support

Indications

- Acute neck injuries
- Chronic cervical syndromes

Ordering information

Sizing Chart

Ref. No.	Size	Measurements	Items per Unit
72819-25	X-S	Length 26–33cm, Height 6.5cm	bx/1
72819-26	S	Length 30–36cm, Height 7.5cm	bx/1
72819-27	M	Length 33–39cm, Height 9cm	bx/1
72819-28	L	Length 36–42cm, Height 10cm	bx/1
72819-29	XL	Length 39–45cm, Height 11cm	bx/1
72819-30	XL (short)	Length 38–48cm, Height 7cm	bx/1

Actimove® Vertebrace

Extrication Collar



- Firm cervical support
- Easy to apply
- Convenient colour coded identification system

Actimove® Vertebrace is a simple to use, lightweight, padded support device used to help immobilise the cervical vertebrae. Actimove® Vertebrace is designed to aid in supporting and limiting movement of the patient's head and neck during transfer to a backboard. Properly applied, the collar completely encircles the neck and provides rigid support for the head, chin and base of the skull. All materials are radiolucent so that x-rays may be taken with the collar in place.

Special Features:

- Large tracheostomy opening
- X-ray transparent
- One piece construction with integral flip-up chin support. No assembly required
- Reinforced back provides firm cervical support
- Stack flat or roll tight for efficient storage



Easy to use sizing guide

Colour-coded Velcro band closure allows quick identification of a full range of sizes



Indications

- Temporary management of traumatized cervical spine to help prevent further injury
- Osteoarthritis and accompanying trauma or strains
- Muscle spasms and musculoligamentous tears or strains
- Cervical tendinitis
- Post cervical fusion
- Foundation for cervical casting

Ordering information

Ref. No.	Size	Items per Unit
Paediatric Short - Light Blue		
79300-00	6.5cm x 43cm	1
Paediatric Red		
79305-00	9cm x 45cm	1
Extra Short – Green		
79301-00	7cm x 57cm	1
Short – Yellow		
79302-00	9cm x 57cm	1
Regular – Orange		
79303-00	11.5cm x 57cm	1
Tall – Navy		
79304-00	14cm x 57cm	1

